


This page will be replaced by "FINESIA BL Implant (Prosthesis) System Manual" : page 47.

Variations in temporary abutment size

(Unit: mm)

Product drawing	Platform	Diameter of the compatible implant (Φ)	Diameter (W)	Post height (H)	Cuff height (GH)
 ST R	NP	3.2/3.4	3.9	9.0	0.9
			3.5		2.0
					3.0
					4.0
					5.0
	RP	3.7/4.2	4.1	9.0	0.7
			4.5		2.0
					3.0
					4.0
					5.0
	WP	4.7/5.2	5.2	9.0	0.7
			5.5		2.0
3.0					
4.0					
5.0					

1-1 Precautions for the use of temporary abutments

- 1) As a standard, an abutment is installed at least 3 months after implant placement in the mandible and at least 6 months after implant placement in the maxilla. Before abutment placement, please ensure that the soft tissues and bone have healed.
- 2) This product is not sterilized.
- 3) Obtain X-ray images as needed and confirm that there is no gap between the abutment and implant body.
- 4) Before attaching the abutment, thoroughly clean the inside of the implant body and eliminate foreign material, including blood.
- 5) A lab screw is used for laboratory work. Because the shape of a lab screw is the same as that of an abutment screw, carefully check the product label before use to avoid confusion.
- 6) The tightening torque for a tightened screw is slightly decreased because of initial plastic deformation of the screw itself. Retightening restores the fastening force, brings the components closer, and prevents loosening. Retighten the abutment screw two or three times.
- 7) Use an alternatively available abutment screw at the time of replacement during maintenance.
- 8) Do not store the products in an area exposed to high temperature, humidity, or direct sunlight. To prevent the risk of infection, never use implant bodies dropped on the floor or contaminated by foreign material such as saliva. Please discard such implants.
- 9) When tightening the abutment or operating a dedicated tool, be careful so that the patient does not accidentally swallow them.
- 10) Do not subject the abutment to excessive force, which can deform or damage it.
- 11) Check the package before opening it. Do not use the abutment, if the package is damaged.
- 12) Do not use any tool other than dedicated ones, otherwise you can damage the implant body.
- 13) Assume that the maximum usable time for temporary abutments is 2 months.
- 14) No occlusal force should be applied to the temporary abutment, considering it is a temporary component.
- 15) Read the package inserts before use and thoroughly understand the function, method of use, and surgical procedure.